Malaria and Malaria PLUS DNA Amplification Assays

DNA Amplification Assays for the Detection of *Plasmodium sp.*

REF 280925, 281125

IVD In vitro diagnostic medical device

INTENDED USE

The *illumigene* Malaria and *illumigene* Malaria PLUS DNA amplification assays, performed on the *illumipro-10™*, are qualitative in vitro diagnostic tests for the direct detection of *Plasmodium sp.* DNA in human venous EDTA whole blood specimens from individuals with signs and symptoms of malarial infection. Results from *illumigene* Malaria assays are intended to be used as an aid in the diagnosis of human malaria infection.

illumigene Malaria assays utilize loop-mediated isothermal DNA amplification (LAMP) technology to detect *Plasmodium sp.* DNA by targeting segments of the *Plasmodium* genome. *illumigene* Malaria assays do not distinguish between *Plasmodium* species.

illumigene Malaria and *illumigene* Malaria PLUS are intended for use in hospital, reference or state laboratory settings. The devices are not intended for nonlaboratory point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

The *illumigene* Malaria and *illumigene* Malaria PLUS molecular assays are based on loop-mediated amplification (LAMP)^{1,2} technology. The assays target a region of the *Plasmodium* genome that is conserved across *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae*, and *Plasmodium knowlesi*. The *illumigene* Malaria assays target a 214 base pair (bp) sequence of a *Plasmodium sp.* mitochondrial DNA noncoding region.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian *illumipro-10* Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The *illumigene* Malaria kit includes *illumigene* Malaria Test Devices, Sample Preparation Apparatus IV (SMP PREP IV), *illumigene* Buffer I, and tubes. Whole blood specimens are first treated with Buffer I; cells are disrupted and nucleic acids are released. The lysate is added to SMP PREP IV and filtered into a clean Tube I by gently squeezing. Sample effluent collected after filtration contains nucleic acids for use with the *illumigene* Malaria Test Device.

The *illumigene* Malaria PLUS kit includes *illumigene* Malaria Test Devices, *M-prep* Columns, *illumigene* Buffer I, *M-prep* Buffer II, *M-prep* Buffer III, and ST Tubes. Whole blood specimens are processed using size exclusion chromatography by gravity flow to separate and purify nucleic acids. *M-prep* Columns provide the matrix for purification of nucleic acids while the included buffers facilitate cell lysis and sample elution. Whole blood specimens are first treated with Buffer I; cells are disrupted and nucleic acids are released. Sample material is applied to the *M-prep* Column and allowed to move through the column by gravity, with the largest molecules moving through the column matrix first, or fastest. *M-prep* Buffer II is applied to the *M-prep* Column, facilitating movement of the sample through the column matrix. *M-prep* Buffer II contains a non-reactive colorant that provides a visual indicator of sample progression. *M-prep* Buffer III is applied to the *M-prep* Column after Buffer II has entered the gel matrix. Sample effluent collected after the addition of *M-prep* Buffer III contains nucleic acids for use with the *illumigene* Malaria Test Device. Neither *illumigene* Malaria nor *illumigene* Malaria PLUS require specialized laboratory equipment.

The *illumigene* Malaria Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with *Plasmodium sp.*-specific primers and a CONTROL chamber with human mitochondrial DNA-specific primers. Human mitochondrial DNA in whole blood samples and the human mitochondrial DNA-specific primers in the Test Device CONTROL chambers function as the Internal Control for the assay. During specimen preparation, human mitochondrial DNA is extracted with the *Plasmodium sp.* mitochondrial DNA to allow for parallel processing of target DNA and Control DNA through amplification and detection. The Internal Control monitors DNA extraction, amplification inhibition, assay reagent performance and sample processing effectiveness. The Control target must be amplified and detected in the final reaction or the test is considered invalid and results are not reported.

The *illumipro-10* monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal_{initial}, S) and at the assay Run End (Signal_{final}, S_r). The *illumipro-10* calculates the change in light transmission between Run End and Run Start (S_r:S) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S_r:S ratios less than 70% are reported as 'POSITIVE'; TEST chamber S_r:S ratios greater than or equal to 70% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_r:S ratios less than 85% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S_r:S ratios greater than or equal to 85% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

Malaria is a public health challenge world-wide, with an estimated 3.2 billion people at risk for infection in 97 countries.³ The World Health Organization (WHO) estimates 198 million cases occurred in 2013 that led to approximately 584,000 deaths.³ Malaria is caused by the intracellular parasite *Plasmodium* and is transmitted to humans by female *Anopheles* mosquitoes.³ Once humans are infected, *Plasmodium* exhibit multiple life stage morphologies for initial growth in liver cells, followed by invasion of red blood cells.⁴ *Plasmodium* replicate within red blood cells, causing cell rupture and clinical symptoms of the disease.⁴

There are five *Plasmodium* species known to cause malaria in humans (*Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium malariae*, *Plasmodium ovale*, and *Plasmodium knowlesi*), with the majority of infections caused by *P. falciparum* and *P. vivax*. Malaria treatment is dependent on the *Plasmodium* species, clinical status of the patient, and antimicrobial susceptibility.⁵ Infections caused by *P. falciparum* and *P. knowlesi* are typically more severe and require rapid treatment.⁵

The WHO recommends diagnostic testing for all individuals with suspected malaria.³ Thin and thick smear microscopy is currently the gold standard diagnostic method. However, accuracy and precision of microscopy is highly dependent on the technician skill (limit of detection ranging from 4-100 parasites/ μ L or higher), with false negative and false positive rates between 10-25% and approximately 2-3 fold variability in quantitation.^{6,7,8} In addition, errors in the *Plasmodium* species identification are common in routine microscopy examination.⁸ Immunooassay-based rapid diagnostic tests (RDTs) are also an available diagnostic method. Quality assessment testing performed by the WHO demonstrated that at low parasite densities (200 parasites/ μ L), only 76% and 42% of available RDTs could detect *P. falciparum* and *P. vivax*, respectively.³

The *illumigene* Malaria assays provide rapid results with higher sensitivity than microscopy and RDTs.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

For *illumigene* Malaria, Catalog Number 280925

- illumigene* Malaria Test Device: Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either *Plasmodium sp.*-specific primers (TEST Chamber) or human mitochondrial DNA-specific primers (CONTROL Chamber).
- illumigene* Buffer I: Lysis solution containing 0.2N sodium hydroxide.
- illumigene* Sample Preparation Apparatus IV (SMP PREP IV): Tris-buffer solution containing 0.09% azide as a preservative.
- illumigene* Tube I: 1.5 mL tubes

For *illumigene* Malaria PLUS, Catalog Number 281125

- illumigene* Malaria Test Device: Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either *Plasmodium sp.*-specific primers (TEST Chamber) or human mitochondrial DNA-specific primers (CONTROL Chamber).
- illumigene* Buffer I: Lysis solution containing 0.2N sodium hydroxide.
- M-prep* Buffer II: Solution containing phenol red and 0.09% azide as a preservative.
- M-prep* Buffer III: Tris-buffered solution with 0.09% azide as a preservative.
- M-prep* Column: 5.0 cm long column with twist-off tip and plug seal cap containing chromatography resin in a tris-buffered solution containing sodium azide (0.09%) as a preservative.
- ST Tubes: 2.0 mL screw-top microcentrifuge tubes.

MATERIALS PROVIDED SEPARATELY

- illumigene* Malaria External Control Kit, Meridian Bioscience, Inc. Catalog Number: 279970

MATERIALS NOT PROVIDED

- Disposable latex gloves, powder free
- DNase/RNase-free, aerosol resistant pipette tips
- Blood collection tubes with EDTA anticoagulant

EQUIPMENT NOT PROVIDED

- Interval Timer
- Vortex mixer (optional)
- Micropipette capable of dispensing 50 μ L
- Micropipette capable of dispensing 250 μ L (Catalog Number 281125 only)
- illumipro-10™*, Meridian Bioscience, Inc. Catalog Number: 610172

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Do not interchange kit reagents and Test Devices between lots. Individual lots of Tube I and ST Tubes are interchangeable within each kit type provided they are within the assigned expiration date when used.
- The *illumipro-10* may produce incorrect results if the Malaria assay program is not used for testing.
- Follow Biosafety Level 2 and Good Laboratory practices during testing.⁹ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
- Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
- Quality Control Programs for Molecular Testing Laboratories, including proper use and care of equipment, should be employed.¹⁰
- The *illumigene* Malaria Test Devices contain lyophilized reagents. The protective pouches should not be opened until ready to perform the assay.
- The *illumigene* Malaria Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
- Dispose of used *illumigene* Malaria Test Devices, *M-prep* Columns, and tubes immediately after processing. Leave the Test Device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD AND PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com, for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit components at the temperature indicated on the label.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (19-30 C) before use. Ensure Buffer I has been warmed to room temperature completely prior to use and no precipitate is visible. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

COLUMN PREPARATION (For *illumigene* Malaria PLUS Catalog Number 281125 only)

Visually inspect each *M-prep* Column and ensure the column bed and filters were not disturbed during transportation. Columns with disturbed resin beds or improperly placed filters may not perform correctly.

Column Set-up:

- Use 1 *M-prep* Column for each sample to be tested. Remove the top cap and the bottom twist-off tip.
- Place the column tip into an ST Tube. The ST Tube should be held at a slight angle (approximately 15 degrees) from the *M-prep* Column.
NOTE: Do not seat the column directly in the collection tube as this may create a seal and cause improper flow through the column. Refer to Figure 1 for proper column and collection tube placement.
- Let the *M-prep* Column drain into the ST Tube. Drained *M-prep* Columns should be used within 1 hour.
- Proceed with Sample Preparation for *illumigene* Malaria PLUS.

FIGURE 1: *M-prep* Column and ST Tube placement

SPECIMEN COLLECTION AND PREPARATION

Sample type: Human venous whole blood samples with EDTA as a preservative.

Sample Collection: Venous whole blood samples should be collected into a specimen tube containing EDTA according with institutional guidelines for collection of clinical specimens for malaria infection. Ensure the blood collection tube is inverted immediately after collection at least 8 times, or according to the manufacturer's instructions.

EDTA venous whole blood samples may be stored at 2-30 C after collection and during transportation to the laboratory. Samples should be tested as soon as possible, but may be stored for up to 7 days at room temperature (19-30 C) or up to 14 days refrigerated (2-8 C) prior to testing. Samples that will not be tested within this time frame should be frozen immediately at -20 C for up to 30 days until tested. Samples may be frozen and thawed 2 times after storage at -20 C prior to testing with the *illumigene* Malaria assays.

SPECIMEN PREPARATION

NOTE: Ensure that the *illumipro-10* is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the *illumipro-10* Operator's Manual for further information regarding instrument set-up and operation.

NOTE: Ensure specimens are at room temperature (19-30 C) before specimen preparation.

Specimens to be tested with the *illumigene* Malaria assay may be prepared using one of the following methods:

1. *illumigene* Malaria Specimen Preparation (Catalog Number 280925)

- Invert EDTA whole blood specimen 2-3 times to mix.
- Add 50 µL of the collected venous whole blood sample (with EDTA) to one tube of Buffer I. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Hold the sample for 2 minutes.
- Mix by inversion 5 times or by vortexing for approximately 10 seconds and immediately transfer 50 µL of lysate into SMP PREP IV. Mix by inversion 5 times or by vortexing for 10 seconds.
- Gently squeeze the SMP PREP IV and slowly collect 5 to 10 drops into a clean Tube I. Visually verify that the eluate is tinted red to reddish brown. Label the tube with the specimen identification and proceed to the Test Procedure.

2. *illumigene* Malaria PLUS Specimen Preparation (Catalog Number 281125)

NOTE: Sample elution steps with *M-prep* Columns should take no longer than 30 minutes. Samples that take longer than 30 minutes to elute should be discarded and re-tested with the original patient sample.

- Invert EDTA whole blood specimen 2-3 times to mix.
- Add 50 µL of the collected venous whole blood sample (with EDTA) to one tube of Buffer I. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Hold the sample for 2 minutes.
- Mix by inversion 5 times or by vortexing for approximately 10 seconds and, using a micropipette, immediately transfer 250 µL of the prepared sample to the top of an appropriately labeled and prepared *M-prep* Column. Wait approximately 2 minutes, or until the sample has been absorbed by the column and flow stops. This step should take no longer than 30 minutes.
- Using a micropipette, add 250 µL of *M-prep* Buffer II to the top of the *M-prep* Column. Discard the pipette tip. The column will have a red appearance after the addition of *M-prep* Buffer II. Wait approximately 2 minutes, or until the red-colored buffer is absorbed by the column and flow stops. This step should take no longer than 30 minutes.
- Remove the last drop of liquid from the column tip with the ST Tube. Discard the tube.
- Place a clean ST Tube under the *M-prep* Column. Using a micropipette, add 250 µL of *M-prep* Buffer III to the top of the *M-prep* Column. Discard the pipette tip. Wait approximately 2 minutes or until flow stops. This step should take no longer than 30 minutes.
- Remove the last drop of liquid from the column tip with the ST Tube. Visually verify that the eluate is tinted red to reddish brown. Label the tube with sample identification information and proceed to the Test Procedure.

TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single *illumipro-10* run.

- Remove 1 *illumigene* Malaria Test Device from its protective pouch per sample. Carefully open the device, holding the chambers such that the lyophilized reagents will not fall out upon opening. Place the device on a flat surface or in a rack that can accommodate the device.
- Using a micropipette, transfer 50 µL of the sample to both the TEST (White Bead) and CONTROL (Yellow Bead) chambers of the *illumigene* Malaria Test. Take care not to introduce air to the reaction mixture.
- Close the *illumigene* Test Device and fasten the latches securely.
- Tap device(s) on the bench top or mix to remove air bubbles. Carefully examine the Test Device(s) for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
- Repeat Test Procedure Steps for all samples to be tested.
- Insert the *illumigene* Test Devices into the *illumipro-10* and initiate run using the Malaria Program. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains <i>Plasmodium</i> sp. target DNA.
	NEGATIVE	No <i>Plasmodium</i> sp. DNA detected.
	INVALID	No reportable result. Repeat the test using the original sample. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No <i>illumigene</i> Test Device in the <i>illumipro-10</i> Well. OR The <i>illumigene</i> Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- Each device contains an internal control that controls for amplification inhibition, assay reagents, DNA preparation, and sample processing effectiveness. Human mitochondrial DNA, which serves as the internal control DNA, is isolated from the whole blood sample and processed through all steps of the procedure. Primers for amplification of internal control DNA are present in the Control Chamber of the *illumigene* Test Device.
- Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
- illumigene* Malaria External Positive Control Reagent is supplied separately (Catalog 279970). Alternatively, previously characterized clinical or contrived *Plasmodium* sp. positive blood samples can be used as an external positive control. A qualified negative human whole blood sample with EDTA may be used as an external negative control. It is recommended that the reactivity of each new lot and each new shipment of *illumigene* Malaria and *illumigene* Malaria PLUS be verified on receipt and before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The *illumigene* Malaria and *illumigene* Malaria PLUS test kit should not be used in patient testing if the external controls do not produce the correct results.
- A separate Test Device must be used for each external control reagent.

EXPECTED VALUES

The incidence of *Plasmodium* sp. as detected by the *illumigene* Malaria assay during the 2015 clinical study was 68.1% (147/216) using the *illumigene* Malaria sample preparation method (Catalog number 280925) and 70.6% (149/211) using the *illumigene* Malaria PLUS Method (Catalog number 281125). The overall incidence of each *Plasmodium* species is provided below.

Species	Prevalence by <i>Plasmodium</i> species		<i>illumigene</i> Malaria PLUS (n=211)	
	Total Positive	Prevalence	Total Positive	Prevalence
<i>P. falciparum</i>	137	63.4%	134	63.5%
<i>P. vivax</i>	0	0%	0	0%
<i>P. ovale</i>	1	0.5%	1	0.5%
<i>P. malariae</i>	1	0.5%	1	0.5%
<i>P. knowlesi</i>	0	0%	0	0%
Unknown*	8	3.7%	13	6.2%

*Samples were negative by microscopy and therefore species could not be determined.

LIMITATIONS OF THE PROCEDURE

- This product can only be used with the *illumipro-10* instrument.
- illumigene* Malaria and *illumigene* Malaria PLUS do not distinguish between *Plasmodium* species.
- Hematocrit (>57.5%) and hemoglobin (>30 g/dL) above normal physiological levels may produce invalid results with the *illumigene* Malaria assays.
- Performance of the *illumigene* Malaria assays has not been established for packed red blood cell specimens.
- Performance of *Plasmodium knowlesi* with the *illumigene* Malaria assay was established using purified genomic DNA only; whole organism testing was not performed.
- illumigene* Malaria and *illumigene* Malaria PLUS are qualitative assays and do not provide quantitative values or information about organism load.
- The detection of nucleic acids is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
- Organism nucleic acid may persist *in vivo*, independent of organism viability. *illumigene* Malaria and *illumigene* Malaria PLUS do not distinguish between viable and nonviable organisms.
- As with all molecular based diagnostic tests, (A) False negative results may occur from the presence of inhibitors, technical error, sample mix-up or low numbers of organisms in the clinical specimen; (B) False positive results may occur from the presence of cross-contamination by target organisms, their nucleic acids or amplified product, and from non-specific signals.

SPECIFIC PERFORMANCE CHARACTERISTICS

The performance characteristics of the *illumigene* Malaria and *illumigene* Malaria PLUS DNA Amplification Assays were established in clinical studies conducted in November 2015 in Senegal, Africa. Performance characteristics of the assay were compared to the reference method, thin and thick film microscopy. *Plasmodium* species identification for all positive samples was determined by microscopy.

Samples included prospective and retrospective venous whole blood specimens with EDTA from patients with signs and symptoms of *Plasmodium* infection. A total of 216 eligible, de-identified whole blood specimens were evaluated. 66 specimens were tested prospectively and 150 specimens were prospectively-collected and stored frozen prior to testing by the *illumigene* Malaria assay (retrospective). Each specimen was prepared by both *illumigene* Malaria and *illumigene* Malaria PLUS sample preparation methods. An additional five samples were considered ineligible for the *illumigene* Malaria PLUS assay due to process error. Therefore, there are 211 eligible samples for *illumigene* Malaria PLUS. All samples were tested prospectively by microscopy.

Analysis of assay performance with prospective and retrospective specimens indicated there is no difference in performance for fresh and frozen specimens with the *illumigene* Malaria and *illumigene* Malaria PLUS assays. The tables below summarize the performance of the *illumigene* Malaria and *illumigene* Malaria PLUS DNA Amplification Assays for both prospective and retrospective samples combined.

illumigene Malaria Performance Compared to Microscopy

	Microscopy			<i>illumigene</i>	Performance				
	Pos	Neg	Total		INV ^a	Sensitivity	100.0%	139/139	97.3-100.0%
<i>illumigene</i> Malaria	Pos	139	8	147	0(1)				
	Neg	0	67	67	2(2)	Specificity	89.3%	67/75	80.3-94.5%
	Total	139	75	214	2				

^aInitial invalid results are reported within the parentheses. The final number of invalid samples remaining after repeat testing is shown before the parenthesis.

illumigene Malaria PLUS Performance Compared to Microscopy

	Microscopy			<i>illumigene</i>	Performance				
	Pos	Neg	Total		INV	Sensitivity	100.0%	136/136	97.3-100.0%
<i>illumigene</i> Malaria PLUS	Pos	136	13	149	0				
	Neg	0	62	62	0	Specificity	82.7%	62/75	72.6-89.6%
	Total	136	75	211	0				

Samples were collected from patients ranging in age from 4 to 77 years. 15 were between 1 – 12 years of age; 75 were between 13 – 21 years of age; 125 were ≥22 years old; the age of one patient was not defined. There was no noted performance difference based on age. The study population included 84 (38.9%) female and 131 (60.6%) male patients; the gender of one patient was not defined. There is no expectation that assay performance is influenced by gender.

ANALYTICAL SENSITIVITY

Limit of Detection (LoD) was determined using *Plasmodium falciparum* and *Plasmodium vivax*. The LoD was confirmed using 20 replicates and a stated probability (eg, 95%, where 19/20 replicates are positive) of obtaining positive responses.

	<i>illumigene</i> Malaria	<i>illumigene</i> Malaria PLUS
<i>Plasmodium</i> Species	Parasites/µL	Parasites/µL
<i>P. falciparum</i> (3D7)	2	0.25
<i>P. vivax</i> (India VII)	0.125	0.063

ASSAY REACTIVITY

Quantified stocks of *P. malariae*, *P. ovale*, and *P. knowlesi* DNA were diluted in negative whole blood to approximately three times the limit of detect for each *illumigene* Malaria and *illumigene* Malaria PLUS sample preparation method (approximately 72 copies/µL or 6 parasites/µL for *illumigene* Malaria; approximately 9 copies/µL or 0.75 parasites/µL for *illumigene* Malaria PLUS). All species reacted at the concentrations tested.

REPRODUCIBILITY

Reproducibility studies were carried out by three internal laboratories. Blind-coded panels of 10 samples were supplied to participating laboratories. The panels included contrived *Plasmodium falciparum* (strain 3D7) samples manufactured as moderate positive samples, low positive samples, or high negative samples. The panel also included one negative whole blood sample. Testing was performed by at least 2 different operators at each laboratory on the same day (intra-assay variability) for five days (inter-assay variability). Three lots each of *illumigene* Malaria or *illumigene* Malaria PLUS and 8 *illumipro-10* instruments were used in this study. External Positive and Negative Controls were tested with each panel; a qualified negative blood sample was used as the Negative Control.

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Reproducibility Study Summary: <i>illumigene</i> Malaria Sample Preparation Method								
Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
Moderate Positive (8 parasites/μL)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive (2 parasites/μL)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
High Negative (0.458 parasites/μL)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Negative	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

Reproducibility Study Summary: <i>illumigene</i> Malaria PLUS Sample Preparation Method								
Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
Moderate Positive (1 parasite/μL)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive (0.5 parasites/μL)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
High Negative (0.029 parasites/μL)	29/30	96.7%	28/30	93.3%	30/30	100.0%	87/90	96.7%
Negative	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

CROSSREACTIVITY STUDIES

Crossreactivity studies employed positive (*P. falciparum* 3D7) and negative whole blood specimens inoculated with bacterial or fungal organisms to a minimum concentration of 1.0×10^6 CFU/mL, virus at a minimum of 1.0×10^5 TCID₅₀/mL, or protozoans to a minimum concentration of 1.0×10^5 organisms/mL. Where whole organisms were not available, 1.0×10^6 copies/mL for genomic DNA was tested. None of the following organisms or their genetic material reacted with the *illumigene* Malaria assays:

Babesia microti, *Borrelia burgdorferi*, *Haemophilus influenzae*, *Klebsiella pneumoniae*, *Leishmania donovani*, *Leishmania infantum*, *Leptospira interrogans*, *Mycobacterium smegmatis*, *Salmonella typhi*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Streptococcus pneumoniae*, *Treponema pallidum*, *Trypanosoma cruzi*, *Trypanosoma rangeli*, *Vibrio parahaemolyticus*, *Cytomegalovirus* (CMV), Epstein-Barr virus (EBV), Hepatitis B virus, Hepatitis C virus, Herpes simplex virus 1 (HSV 1), HIV-1, Human Papilloma Virus (HPV), and Rubella virus.

The following organisms or their DNA could not be obtained for testing and were evaluated through *in silico* analysis. None of the following organisms are expected to react with the *illumigene* Malaria assays based on this analysis: *Anaplasma phagocytophilum*, *Clostridium botulinum*, *Orientia tsutsugamushi*, *Rickettsia conorii*, *Rickettsia rickettsii*, *Rickettsia typhi*, Chikungunya virus, Dengue virus (types 1-4), West Nile virus, and Yellow Fever virus.

Human genomic DNA was nonreactive at 1.0×10^6 copies/mL.

TESTS FOR INTERFERING SUBSTANCES

Interference testing was performed in the presence of worst-case concentrations of chemical and biological substances introduced directly into contrived low positive (*P. falciparum* 3D7) and negative whole blood samples. No interference was observed with the following substances for either *illumigene* Malaria sample preparation methods:

Acetaminophen (1 mg/mL), Amoxicillin (0.1 mg/mL), Artesunate (1 mg/mL), Aspirin (1 mg/mL), Atovaquone (0.1 mg/mL), Cephalaxin (0.1 mg/mL), Chloroquine (1 mg/mL), Ciprofloxacin (0.1 mg/mL), Clindamycin (1 mg/mL), Doxycycline hydrochloride (1 mg/mL), Erythromycin (0.1 mg/mL), Hydroxychloroquine sulfate (1 mg/mL), Ibuprofen (1 mg/mL), Lumefantrine (1 mg/mL), Mefloquine (1 mg/mL), Primaquine phosphate (1 mg/mL), Proguanil hydrochloride (1 mg/mL), Pyrimethamine (1 mg/mL), Quinine sulfate (1 mg/mL), Sodium Citrate (0.11 M), elevated Bilirubin (>0.15 mg/dL), elevated leukocytes (buffy coat) >10% v/v, serum albumin (>0.03 g/mL), and triglycerides >9.9 mg/dL.

Hematocrit (>57.5%) and hemoglobin (>30 g/dL) above normal physiological levels may produce invalid results with the *illumigene* Malaria assays.

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REFERENCES

- Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. *Mol Cell Probes* 2002;16:223-29.
- Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. *J Biochem Biophys* 2004;59:145-47.
- World Malaria Report 2014. Global Malaria Programme, World Health Organization. http://www.who.int/malaria/publications/world_malaria_report/en/.
- Cowman AF and Crabb BS. Invasion of Red Blood Cells by Malaria Parasites. *Cell* 2006; 124: 755-766.
- Treatment of Malaria (Guidelines for Physicians). Centers for Disease Control and Prevention, July 2013. www.cdc.gov/malaria/resources/pdf/clinicalguidance.pdf.
- O'Meara WP, Barcus M, Wongsrichanalai C, Muth S, Maguire JD, Jordan RG, Prescott WR, and McKenzie FE. Reader Technique as a Source of Variability in Determining Malaria Parasite Density by Microscopy. *Malaria J* 2006;5:118.
- McKenzie FE, Sirichaisinthop J, Miller RS, Gasser RA, and Wongsrichanalai C. Dependence of Malaria Detection and Species Diagnosis by Microscopy on Parasite Density. *Am J Trop Med Hyg* 2003;69(4):372-376.
- Wongsrichanalai C, Barcus MJ, Muth S, Sutarnhardja A, and Wernsdorfer WH. A Review of Malaria Diagnostic Tools: Microscopy and Rapid Diagnostic Test (RDT). *Am J Trop Med Hyg* 2007;77(6):119-127.
- US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
- CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute, 2006.

SN1102

REV. 01/16

Meridian Bioscience, Inc.
 USA/Corporate Office
 3471 River Hills Drive
 Cincinnati, Ohio 45244
 Telephone: 513.271.3700
 Orders/Customer Service:
 800.543.1980
 Technical Support Center:
 800.343.3858
 Information Fax: 513.272.5432
 Ordering Fax: 513.271.0124

Meridian Bioscience Europe S. r. l.
 Via dell'Industria, 7
 20020 Villa Cortese, Milano
 ITALY
 Tel: +39 0331 43 36 36
 Fax: +39 0331 43 36 16
 Email: info@meridianbioscience.eu
 WEB: www.meridianbioscience.eu

EC REP
 Authorized Representative

Meridian Bioscience Europe s.a./n.v.
 2 Avenue du Japon - 1420 Braine l'Alleud
 BELGIUM
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 Email: info.bn@meridianbioscience.eu





Meridian Bioscience Europe France
 34 rue de Ponthieu - 75008 Paris
 FRANCE
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
 Postbus 301 - 5460 AH Veghel
 NETHERLANDS
 Tel: +31 (0) 411 62 11 66
 Fax: +31 (0) 411 62 48 41
 Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusqu'à / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
CE	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In-Vitro Diagnostica 98/79/EG	SMP PREP OIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
		STERILE R	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Sterilisation par irradiation aux rayons gamma / Esterilizado por irradiación gamma / Sterilisation durch Gammastrahlen
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	RoHS	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		Caution, consult accompanying documents / Attention, voir les instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	STERILE EO	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Sterilisation par oxyde d'éthylène / Esterilizado por óxido de etileno / Sterilisation durch Ethylenoxid
	Contains sufficient for "n" tests / Contiene sufficiente per "n" test / Contenu suffisant pour "n" tests / Contenido suficiente para "n" ensayos / Inhalt ausreichend für "n" Prüfungen	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction / tamponee / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Límite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform geprüfte
	Female / Femminile / De sexe féminin / Hembra / Frau		Recycle - do not dispose of as general waste / Riciclare - non eliminare come rifiuto generico / Recycler - ne pas jeter dans une poubelle / Recycle - no desecho como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
	Single Use Only / Prodotto Monouso / A usage unique / Para Uso Solo / nur für die einmalige Anwendung		
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	HT TUBE	Heat Treatment Tube / Provetta per il trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhrchen zur Hitzebearbeitung
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Reactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Eviter toute exposition au faisceau / Radiación Laser: Evite la Exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden		HOT SURFACE: Keep hands away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	CAUTION: Laser Radiation / AVERTISSEMENT: Radiation Laser / AVERTISSEMENT: Rayonnement laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung	IPX-0	CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risque de danger / Precaución: Peligroso / WARNUNG: Risikofaktor	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
MIN OIL	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl	MEDIA	Media / Terreno di trasporto / Milieux / Medio / Medium
ST TUBE	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhrchen mit Schnappverschluss	COL	Sample Preparation Column / Colonna di preparazione dei campioni / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufbereitung
BUF SMP	Sample Buffer / Soluzione tampone per il campione / Tampone de l'échantillon / Tampón de muestra / Probenpuffer	PRE REAG	Pretreatment Reagent / Reagente di Pretattamento / Réactif de prétraitement / Reactivo de prettreatment / Reagenz für die Vorbehandlung
		SMP PREP	Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenvorbereitung

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.