

illumipro-10™

Automated isothermal amplification and detection system that provides fast, accurate results in a flexible package that fits into any lab.



- Optimal flexibility for stat or batch runs
- True walk away capability for optimized workflow
- Very low maintenance requirements enhance efficiency
- Instrument verifications can be performed with in-house personnel

illumipro-10™

Digital display guides users through process

Automatic system self-check

Barcode reader for sample ID entry

Reaction blocks operate independently


Small footprint of instrument is only 8.5" x 11"



Procedure



1. Load samples



2. Press "RUN"



3. Enter patient ID

Ordering information

<i>illumipro-10™</i> Incubator Reader.....	610172
<i>illumipro-10™</i> Printer.....	610173
<i>illumipro-10™</i> Keyboard.....	610174
<i>illumipro-10™</i> Operations Manual.....	11007

i For more information, contact an *illumigene*® specialist at 1-888-763-6769 or visit us on the web at www.meridianbioscience.com.